



Disposable Bulldog Clip

Atraumatic traction – No Excessive Occlusion

Advantages:

- Features great traction without increase of occlusion pressure
- Unique jaw cover pad distributes occlusion force over larger contact area ensuring atraumatic traction
- Soft silicone jaw cover pads 3mm wide and 18mm long
- Radiopaque – easy identification when needed

Materials:

- Polycarbonate body
- Medical grade stainless steel spring
- Soft, silicone jaw cover pads
- Radiopaque

Jaw Design:

- Straight, parallel
- Low profile
- Jaw covers – 3mm wide & 18mm long
- Max. opening 6mm



CV 700



Cardio Vision[®]

MIC-Disposable Bulldog Clip

| Name | Description | Item Number |
|----------------------------|---|-------------|
| CV Disposable Bulldog Clip | soft silicone jaw cover pads 3mm wide, 18mm long | CV700 |